

Assay Kits

Heat Shock


StressMarq
 Biosciences INC

Antibodies • Proteins • Kits • Small Molecules

Heat Shock StressXpress® ELISA Kits

StressMarq offers a range of Enzyme-Linked Immunosorbent Assay (ELISA) kits in the area of cellular stress, specifically for Heat Shock Protein (HSP) research. Our kits are rigorously tested to ensure they meet our high standards for quality and reproducibility. All kits are available in two sizes: the standard 96 well size and 5 x 96 well size for high throughput screening.

Product Name	Species Reactivity	Sample Type	Sensitivity	Catalog No.
Alpha B-Crystallin StressXpress® ELISA Kit	Human	Cell lysates Tissue extracts Serum	0.009 ng/mL	SKT-123 SKT-123-480
HSP27 StressXpress® ELISA Kit	Human	Cell lysates Tissue extracts Serum	0.04 ng/mL	SKT-109 SKT-109-480
HSP60 StressXpress® ELISA Kit	Human	Cell lysates Tissue extracts Serum	1.37 ng/mL	SKT-110 SKT-110-480
HSP65 StressXpress® ELISA Kit	M. bovis M. tuberculosis	Cell lysates Tissue extracts Serum	0.06 ng/mL	SKT-113 SKT-113-480
HSC70 StressXpress® ELISA Kit	Human	Cell lysates Tissue extracts Serum	1.54 ng/mL	SKT-106 SKT-106-480
HSP70 StressXpress® ELISA Kit	Human	Cell lysates Tissue extracts	0.18 ng/mL	SKT-105 SKT-105-480
High-Sensitivity HSP70 StressXpress® ELISA Kit	Human Monkey Canine Rat & Mouse	Serum Plasma Cell lysates Tissue extracts	0.02 ng/mL	SKT-108 SKT-108-480
Plant HSP70 StressXpress® ELISA Kit	Flax Alfalfa Corn Human	Plant extract Cell lysates Tissue extract Serum & Plasma	0.18 ng/mL	SKT-124 SKT-124-480
HSP90α StressXpress® ELISA Kit	Human	Cell lysates Tissue extracts Serum	0.5-32ng/mL	SKT-107 SKT-107-480