

# Oxidative Stress Kits



Antibodies • Proteins • Kits • Small Molecules

## Oxidative Stress StressXpress® Assay Kits

StressMarq offers a range of ELISA kits and Activity Assay kits in the area of oxidative stress. Our kits are rigorously tested to ensure they meet our high standards for quality and reproducibility.

Product Name	Species Reactivity	Sample Type	Sensitivity	Catalog No.
<b>Catalase Activity</b> StressXpress® Colorimetric Assay Kit	Species Independent	Serum, Plasma, Cells, Tissues, Erythrocyte Lysates	0.052 U/mL	SKT-215
<b>DNA Damage</b> StressXpress® EIA kit	Species Independent	Urine Cell culture Plasma Sample matrices	30pg/mL	SKT-120
<b>Glutathione</b> StressXpress® Fluorescent Assay Kit	Species Independent	Whole Blood, Serum, Plasma, Erythrocytes, Urine, Cell Lysates, Tissue Samples	45 nM Free 48 nM Total	SKT-202
<b>Glutathione Reductase Activity</b> StressXpress® Fluorescent Assay Kit	Species Independent	Lysate RBCs Serum Plasma	0.009 mU/mL	SKT-204
<b>Glutathione S-Transferase Activity</b> StressXpress® Fluorescent Assay Kit	Species Independent	Lysates Serum Plasma	2.70 mU/mL	SKT-203
<b>HO-1</b> StressXpress® ELISA Kit	Human	Cell lysates Tissue extracts Serum	0.21ng/mL	SKT-111
<b>Hydrogen Peroxide</b> StressXpress® Colorimetric Assay Kit	Species Independent	Fresh Urine Buffer TCM	1.83 µM	SKT-216
<b>SOD Activity</b> StressXpress® Colorimetric Assay Kit	Species Independent	Serum Urine Buffer	0.044 U/mL	SKT-214